

# Disposable Haemorrhoid Stapler

Sterile Disposable Stapler has been designed for the surgical procedure for prolapse and haemorrhoids (PPH) throughout the anal canal.

Allows controlled tissue compression and adjustable-height staples for uniform staple formation at the specific height. The anvil is detachable for better visual control.

Clear Anoscope offer reference markings for enhanced visualisation of the anatomy and precise placement of the purse string sutures.

Controlled tissue compression and adjustable-height staples for uniform staple formation.

Closed staple height of 0.75mm-1.5mm – designed to increase staple line compression on tissue and key blood vessels.

## Disposable Haemorrhoid Stapler

reference code	Number of Staples	Diameter	Height
WM-CSH-32	30	32mm	4.2mm
WM-CSH-34	32	34mm	4.2mm

