

Disposable Linear Stapler

One-handed design gives improved ergonomics and ease of use.

Disposable Reloadable Linear Stapler has application throughout the alimentary tract and in thoracic surgery for transection and resection of internal tissues.

Available with two staple heights, which facilitate anastomosing of thick and standard tissues.

Disposable Linear Stapler

reference code	cartridge code	closed staple height	cartridge colour
WM-LS-30-38	WM-LSAC-30-38	1.5mm	Blue
WM-LS-30-48	WM-LSAC-30-48	2.0mm	Green
WM-LS-45-38	WM-LSAC-45-38	1.5mm	Blue
WM-LS-45-48	WM-LSAC-45-48	2.0mm	Green
WM-LS-60-38	WM-LSAC-60-38	1.5mm	Blue
WM-LS-60-48	WM-LSAC-60-48	2.0mm	Green
WM-LS-90-42	WM-LSAC-90-42	1.5mm	Blue
WM-LS-90-48	WM-LSAC-90-48	2.0mm	Green

