


Polymer Ligating Clips and Applicators



Clips are made of non-absorbable polymer, inert, non-toxic and well tolerated by the body, highly resistant in time.

Flexible design of the arms to allow optimal fit on vessels.

Flexible hinge that allows easy opening and closing.

Integrated teeth to ensure greater stability and prevent slippage on the tissue.

Distinct tactile feedback confirms jaw seating and secure vessel placement.

Bosses are designed to securely retain clip in applicator jaws.

Applicators for open and endoscopic surgery equipped with precision anchoring system that allows clips to remain stable without risk of accidental loss.

Loaders and applicators with code colour for easy identification.

Non-Absorbable polymer is inert, nonconductive, radiolucent and does not interfere with CT, MRI, or X-ray diagnostics.

Clips contained in colour coded cartridges for ease of identification.



Our broad range of applier lengths including 45 cm endoscopic appliers for extended reach allows the surgeon to deliver clips with control and confidence

Medium-Large appliers available in 5 mm diameter.

20° angled appliers for effective control and better access.



Take-Apart
2-piece
Endoscopic
Appliers

Polymer Ligating Clips

reference code	size	vessel range	colour	pcs /cartridge	cartridges /box
WM-PC230	Med-Large	3-10mm	Green	6	14
WM-PC240	Large	5-13mm	Purple	6	14
WM-PC250	Extra Large	7-16mm	Yellow	6	14

Endoscopic Appliers

reference code	size	colour code	diameter	length
WM-PC205ML	Medium-Large	Green	5mm	33cm
WM-PC205L	Large	Purple	10mm	33cm
WM-PC205EL	Extra Large	Yellow	10mm	33cm
45cms Extra Length Endoscopic Appliers				
WM-PC205ML45	Medium-Large	Green	5mm	45cm
WM-PC205L45	Large	Purple	10mm	45cm
WM-PC205EL45	Extra Large	Yellow	10mm	45cm
20° Angled Tip Appliers				
WM-PC205MLArt	Medium-Large	Green	5mm	33cm
WM-PC205LArt	Large	Purple	10mm	33cm
WM-PC205ELArt	Extra Large	Yellow	10mm	33cm
Endoscopic Remover				
WM-PCR	Med-Large/Large	Black	5mm	33cm

Endoscopic Appliers 2-Piece

reference code	size	colour code	diameter	length
WM-PC205ML2P	Medium-Large	Green	5mm	33cm
WM-PC205L2P	Large	Purple	10mm	33cm
WM-PC205EL2P	Extra Large	Yellow	10mm	33cm

Open Appliers

reference code	size	colour code	length
Open Appliers			
WM-PCML25	Medium-Large	Green	25cm
WM-PCL25	Large	Purple	25cm
WM-PEL25	Extra Large	Yellow	25cm
Open Appliers 70° Angle			
WM-PCML25A	Medium-Large	Green	25cm
WM-PCL25A	Large	Purple	25cm
WM-PEL25A	Extra Large	Yellow	25cm

20° angled appliers for effective control and better access.

